

SUTURE MATERIAL BRAND NAME CROSS REFERENCE

SUTURAS Surgeasy	ETHICON	B. BRAUN	D & G (COVIDIEN)	USSC (COVIDIEN)	MYCO MEDICAL	CP MEDICAL	VPL	DEME-TECH
PGA, Braided	VICRYL, Coated (polyglactin 910)	SAFIL	DEXON II	POLYSORB	COATED PGA	VISORB Coated	POLYMEND ES	POLYGLYCOLIC ACID / POLYGLACTIN 910
N/A	VICRYL-RAPID, Coated	SAFIL QUICK	N/A	BIOSYN	N/A	VISORB QUICK	N/A	RAPID ABSORBABLE
PCL Mono (poliglecprone)	MONOCRYL (poliglecprone 25)	MONOSYN	BIOSYN	N/A	N/A	MONOSWIFT	MONOMEND MT	POLYIGLECAPRONE 25
PDN Mono (polydioxanone)	PDS II (polydioxanone)	N/A	MAXON	N/A	N/A	MONO-DOX	MONOMEND MaX	POLYDIOXANONE
GUT Plain	SURGICAL GUT PLAIN	SOFTCAT Plain	GUT Plain	SURGIGUT Plain	GUT - Plain	GUT PLAIN	N/A	CATGUT PLAIN
GUT Chromic	SURGICAL GUT CHROMIC	SOFTCAT Chromic	GUT Chromic	SURGIGUT Chromic	GUT - Chromic	GUT CHROMIC	CHROMEND Chromic Gut	CATGUT CHROMIC
Polyester, Coated Braided	ETHIBOND EXCEL Braided Polyester	SYNTHOFIL	TI-CRON	SURGIDAC	N/A	N/A	N/A	POLYESTER (No Cover)
NYLON, Mono	ETHILON Mono Nylon	DAFILON	DERMALON	MONSOF	BLACK NYLON	MONOMID	NY-ST A	NYLON MONO
POLYESTER, Braided	MERSILENE Polyester	DAGROFIL	DACRON	N/A	POLYESTER- GREEN BRAIDED	POLYBOND	MULTI-ST A	POLYESTER
N/A	NUROLON Braided Nylon	SUPRAMID	SURGILON	BRALON		POLYAMID	N/A	NYLON BRAIDED
SILK, Black Braided	PERMA-HAND SILK	SILKAM	SURGICAL SILK	SOFSILK	SILK – Black Braided	SILK	N/A	SILK
POLYPROPYLENE, Mono	PROLENE Mono Polypropylene	PREMILENE	SURGILENE	SURGIPRO	POLYPROPYLENE	POLYPRO	PRO-ST A	POLYPROPYLENE