

# SUTURAS MEDICAL, LLC

## **SURGEASY PDN (Polydioxanone) SYNTHETIC ABSORBABLE SURGICAL SUTURE**

(THE UNITED STATES PHARMACOPEIA)

USP SIZE	METRIC SIZE (EP)	DIAMETER (mm)
7-0	0.5	0.050 – 0.094
6-0	0.7	0.095 – 0.149
5-0	1	0.150 – 0.199
4-0	1.5	0.200 – 0.249
3-0	2	0.2050– 0.339
2-0	3	0.340– 0.399
0	3.5	0.400 – 0.499
1	4	0.500 – 0.570
2	5	0.571 – 0.610

## **SURGEASY PDN (Polydioxanone) SYNTHETIC ABSORBABLE SURGICAL SUTURE**

(THE UNITED STATES PHARMACOPEIA)

USP SIZE	METRIC SIZE (EP)	MINIMUM AVERAGE VALUE (N)
7-0	0.5	1.4
6-0	0.7	2.5
5-0	1	6.8
4-0	1.5	9.5
3-0	2	17.5
2-0	3	26.8
0	3.5	39.0
1	4	50.8
2	5	63.5

# SUTURAS MEDICAL, LLC

## SURGEASY PCL (Poly (glycolic-co-caprolactone)) SYNTHETIC ABSORBABLE SURGICAL SUTURE

(THE UNITED STATES PHARMACOPEIA)

USP SIZE	METRIC SIZE (EP)	DIAMETER (mm)
7-0	0.5	0.050 – 0.094
6-0	0.7	0.095 – 0.149
5-0	1	0.150 – 0.199
4-0	1.5	0.200 – 0.249
3-0	2	0.2050– 0.339
2-0	3	0.340– 0.399
0	3.5	0.400 – 0.499
1	4	0.500 – 0.570
2	5	0.571 – 0.610

## SURGEASY PCL (Poly (glycolic-co-caprolactone)) SYNTHETIC ABSORBABLE SURGICAL SUTURE

(THE UNITED STATES PHARMACOPEIA)

USP SIZE	METRIC SIZE (EP)	MINIMUM AVERAGE VALUE (N)
7-0	0.5	1.4
6-0	0.7	2.5
5-0	1	6.8
4-0	1.5	9.5
3-0	2	17.5
2-0	3	26.8
0	3.5	39.0
1	4	50.8
2	5	63.5